


# Product Specification Sheet

## OUR Claw



### Product Description:

OUR Claw™ helps the Haz Mat Technician interpret the results from reagent paper strips in the hot zone. These results can then be directly communicated to the entry supervisor. This helps to eliminate, "...the second paper from the right just turned green..."

The "dippable" fingers of OUR Claw™ allow a specific test to be performed on a long neck container without cross contamination to the other reagent paper strips. Metal clips allow both a "buffered" and dry reagent strip to be attached to a finger if desired.

### Product Features and Specifications:

- Metal clips make it easy to attach and detach multiple reagent paper strips in the field
- Appropriately sized to be taken into the hot zone and be "readable" for the haz mat technician
- Uses reagent paper strips to detect chemicals
- Reagent paper strips can be folded over and compared to gradient color charts (pH)
- OUR Claw™ is reusable
- "Dippable" legs allow a sample to be obtained from a long neck container

<b>Product Name</b>	OUR Claw™	<b>Size</b>	4" x 10"
<b>Detection methodology</b>	Fluoride, pH, K, M8	<b>Material</b>	5 mm laminate
<b>Shelf Life (papers)</b>	18 months	<b>OUR Claw™</b>	Reusable

## IN OUR GEAR, LLC

Tools for the first responder...  
by the first responder.

[www.inourgear.net](http://www.inourgear.net)